

HUAKUN (BIPAP/NIV CIRCUIT) with Proximal Line

The expiratory filter provides greater than 99.99% bacterial and viral efficiency to protect clinicians. Filtration offers reassurance and peace of mind to any clinician with concerns about the safety and appropriateness of using NIV during a pandemic.



Description	Quantity	Part Number
Disposable Anesthetic Circuit-Corrugated 1.8m (BIPAP/NIV Circuit) with Proximal Line (Bacteria Filter Optional)	50 Pcs/Box	A10024

Specifications:

- Single Limb – Tubing, 22mm ID
- Disposable circuit is for single use
- Consists of:
 - a) Smooth inner lumen tubing
 - b) Proximal pressure line
 - c) Disposable exhalation port
 - d) Hose clips
 - e) Tubing hanger
 - f) Exhaust vent
- Fit NIV vented mask with male/female 22mm connector
- Does not include natural rubber latex or DEHP
- Two types of tubing:
 - a) 180cm
 - b) 90cm (±10cm)
- Fit with all type ventilators.
- Tubing connector adapter Male 22mm Female 22mm connector
- Optimal humidity and able to control condensation.
- Sterile individual pack
- Registered under Medical Device Act