


Yankauer Suction Handles Adult & Pediatric

Our Yankauer Suction Handles are available with straight and crown tips and with and without vents and with an on/off switch. Each comes with a step connector end to ensure a firm connection to any standard suction tubing. They are also available as standard and double wrapped versions.

- ⊗ Smooth and clean edges for atraumatic suctioning
- ⊗ Stepped connectors
- ⊗ Vented and unvented options
- ⊗ Straight and crown tip options
- ⊗ Double wrapped options
- ⊗ Transparent
- ⊗ Paediatric and adult versions



Yankauer Suction Handles

RE-ORDER DESCRIPTION	OUTER TIP DIAMETER	STERILE /NS	CARTON QTY
F12 Yankauer Handle - Fine Tip Paediatric with Vent	4mm	Sterile	100
F14 Yankauer Handle - Fine Tip Paediatric without Vent	4mm	Sterile	100
F16 Yankauer Handle - Straight Tip with Vent	8mm	Sterile	100
F18 Yankauer Handle - Straight Tip without Vent	8mm	Sterile	100
F20 Yankauer Handle - Crown Tip with Vent	8mm	Sterile	100
F24 Yankauer Handle - Crown Tip without Vent	8mm	Sterile	100