


Disposable Anesthesia Circuit-Coaxial



- Less weight than dual limb circuit
- Low airflow resistance directly bonded connectors to the tube for avoiding leakage risks
- All fittings comply with ISO 5356 standard
- Used with ventilator for delivering anesthesia gas
- Divided into adult, pediatric types

Product name	Type	Length	Product code
Disposable Anesthesia Circuit-Coaxial	adult	1.6M	MZG-04-01-1.6M
		1.8M	MZG-04-01-1.8M
		2.0M	MZG-04-01-2M
		2.4M	MZG-04-01-2.4M
		2.8M	MZG-04-01-2.8M
		3.0M	MZG-04-01-3M
	pediatric	1.6M	MZG-04-02-1.6M
		1.8M	MZG-04-02-1.8M
		2.0M	MZG-04-02-2M
		2.4M	MZG-04-02-2.4M
		2.8M	MZG-04-02-2.8M
		3.0M	MZG-04-02-3M