

Disposable Anesthesia Circuit-Corrugated



- Low airflow resistance directly bonded connectors to the tube for avoiding leakage risks
- Bio-compatible and phthalates free material reduces allergy risk in patients
- Standard 22mm/15mm/10mm connection meet ISO 5356-1
- Used with ventilator for delivering anesthesia gas
- Divided into adult, pediatric and neonatal types

Product name	type	Length	Product code
Disposable Anesthesia Circuit-Corrugated	Adult	1.5M	MZG-01-01-1.5M
		1.6M	MZG-01-01-1.6M
		1.8M	MZG-01-01-1.8M
		2.0M	MZG-01-01-2M
		2.4M	MZG-01-01-2.4M
		2.8M	MZG-01-01-2.8M
		3.0M	MZG-01-01-3M
	Adult,detachable	1.5M	MZG-01-01-1.5M-D
		1.6M	MZG-01-01-1.6M-D
		1.8M	MZG-01-01-1.8M-D
		2.0M	MZG-01-01-2M-D
		2.4M	MZG-01-01-2.4M-D
		2.8M	MZG-01-01-2.8M-D
		3.0M	MZG-01-01-3M-D
		1.5M	MZG-01-02-1.5M

SHAOXING UNDIS MEDICAL TECHNOLOGY CO., LTD.

No.12 Yuexing Technology Park, No.11 Fanrong Road, Binhai New Area, 312366 Shaoxing, Zhejiang Province, PEOPLE'S REPUBLIC OF CHINA

	Pediatric	1.6M	MZG-01-02-1.6M
		1.8M	MZG-01-02-1.8M
		2.0M	MZG-01-02-2M
		2.4M	MZG-01-02-2.4M
		2.8M	MZG-01-02-2.8M
		3.0M	MZG-01-02-3M
	Pediatric,detachable	1.5M	MZG-01-02-1.5M-D
		1.6M	MZG-01-02-1.6M-D
		1.8M	MZG-01-02-1.8M-D
		2.0M	MZG-01-02-2M-D
		2.4M	MZG-01-02-2.4M-D
		2.8M	MZG-01-02-2.8M-D
		3.0M	MZG-01-02-3M-D