



Respiratory Tract Infections

Etiology can be Different

ACRO Respiratory Tract Infection Range

Reliable Rapid Tests for :

✓ Respiratory Syncytial Virus (RSV) ✓ *Adenovirus pneumoniae* ✓ Influenza A+B

ACRO provides the rapid tests with high accuracy:

Analyte	Sensitivity	Specificity	Accuracy
Flu A	99%	98.9%	98.9%
Flu B	97.7%	99%	98.6%
RSV	92.7%	98%	95.6%
Adeno	98.2%	98.2%	98.2%



ACRO Respiratory Tract Infection Range

● Influenza

- Influenza, commonly known as “the flu”, spreads around the world in yearly outbreaks, resulting in about three to five million cases of severe.
- Three of the four types of influenza viruses affect people, Type A, Type B, and Type C. Influenza A, B viruses cause significant morbidity and mortality worldwide. Death occurs mostly in the young, the old, and those with other health problems.
- In the northern and southern parts of the world, outbreaks occur mainly in the winter, while around the Equator, outbreaks may occur at any time of the year.


● Respiratory Syncytial Virus (RSV)

- RSV is a major cause of lower respiratory tract infections and hospital visits during infancy and childhood.
- In temperate climates there is an annual epidemic during the winter months. In tropical climates, infection is most common during the rainy season.

● Adenovirus pneumoniae

- Adenovirus is one of the main pathogens of viral pneumonia.
- Infections in children or young adults usually cause upper respiratory tract diseases. Infections in infants can cause severe acute adenoviral pneumonia.

Convenient Procedure



Ordering Information

Cat. No.	Product Description	Specimen	Format	Pack	CE Status
IAD-502	Adenovirus pneumoniae Antigen Rapid Test	Swab	Cassette	20 T	CE
IIN-501	Influenza A+B Rapid Test	Swab / Nasal Aspirate	Dipstick	25 T	CE
IIN-502	Influenza A+B Rapid Test	Swab / Nasal Aspirate	Cassette	20 T	CE
IFLA-501	Influenza A Rapid Test	Swab / Nasal Aspirate	Dipstick	25 T	CE
IFLA-502	Influenza A Rapid Test	Swab / Nasal Aspirate	Cassette	20 T	CE
IRS-502	Respiratory Syncytial Virus (RSV) Rapid Test	Swab / Nasal Aspirate	Cassette	20 T	CE
IRIC-525	RSV&Influenza A+B Combo Rapid Test Cassette	Swab / Nasal Aspirate	Cassette	10 T	CE
IARC-525	Adenovirus&RSV Combo Rapid Test Cassette	Nasopharyngeal Swab	Cassette	10 T	CE
IARI-535	Adenovirus, RSV and Influenza A+B Combo Rapid Test Cassette	Nasopharyngeal Swab	Cassette	10 T	CE