

BACTERIAL FILTER

● Bacterial Filter with Port for Pediatric or Neonatal



PEDIATRIC



NEONATAL



LATEX
FREE

RFE:ZF-158

Dead Space:20ml

Size:

H:58mm; w. 49mm

Weight:15g

Material:

Housing:

Polypropylene (PP)

Filter: Electrical Cotton

Hydrophobic

Connector:

22M/15F-15M

Tidal volume:50ml-300ml

Gas Sampling Port: Yes

BFE:99.9999% @ 0.5Microns

Sterilization: E.O

VFE:99.999% @ 0.027Microns

ZF-158	MOQ:200pcs	Single pack in PE-Paper bag	200pcs/carton	carton size:50x40x32cm	weight: 5kg
ZF-158	Single pack in PE-Paper bag	50pcs/Per Pack	250pcs/carton	carton size:50x40x47cm	weight:6.2 kg