

COMBINED HME & PLEATED HEPA FILTER

● Combined HME and Pleated HEPA Filter



RFE:ZF-170

Size:

H:88mm; w. 75mm L.50mm

Material:

Housing:PP

HME Part: Paper HME

Filter: Pleated HEPA paper

BFE:99.99999% @ 0.5Microns

VFE:99.99999% @ 0.027Microns

Dead Space:81ml

Weight:50g

Connector:

22M/15F-22F/15M

Tidal volume:300ml-1500ml

Gas Sampling Port: Yes

Sterilization: E.O



ZF-170

MOQ:200pcs

Single pack in PE-Paper bag

150pcs/carton

carton size:50x40x47cm

weight: 8.9kg



**ZHENFU
MEDICAL**