

HMEF FILTER

HME with Bacterial/Viral Filter

● HMEF filter with Port for Adult



RFE:ZF-WKR-001

Dead Space:35ml

Size:

H:74mm; w. 57mm

Weight:30g

Material:

Housing: Polypropylene (PP)

Filter: Electrical Cotton

Hydrophobic

Connector:

22M/15F-22F/15M

Tidal volume:100ml-1500ml

Gas Sampling Port: Yes

BFE:99.9999% @ 0.5Microns

VFE:99.999% @ 0.027Microns

Sterilization: E.O

ZF-WKR-001	MOQ:200pcs	Single pack in PE-Paper bag	200pcs/carton	carton size:50x40x32cm	weight: 7.4kg
ZF-WKR-001	Single pack in PE-Paper bag	50pcs/Per Pack	200pcs/carton	carton size:50x40x47cm	weight: 8.6kg