

HMEF FILTER NEONATAL

● HMEF filter with Port for Neonatal



RFE:ZF-0062

Size:

H:51mm; w. 38mm

Material:

Housing:PP

HME Part: Paper HME

Filter: Electrical Cotton

BFE:99.9999% @ 0.5Microns

VFE:99.999% @ 0.027Microns

Dead Space:8ml

Weight:8g

Connector:

15F-15M/8M

Tidal volume: 50ml-150ml

Gas Sampling Port: Yes

Sterilization: E.O

ZF-0062	MOQ:200pcs	Single pack in PE-Paper bag	200pcs/carton	carton size:50x40x32cm	weight: 3.0kg
ZF-0062	Single pack in PE-Paper bag	50pcs/Per Pack	250pcs/carton	carton size:50x40x47cm	weight:4.3 kg

