

# HMEF FILTER PEDIATRIC

## ● HMEF filter with Port for Pediatric



RFE:ZF-002P

**Dead Space:**15ml

**Size:**

H:67mm; w. 44mm

**Weight:**13.5g

**Connector:**

22M/15F-22F/15M

**Material:**

Housing:

Polypropylene (PP)

Filter: Electrical Cotton

Hydrophobic

**Tidal volume:**150ml-300ml

**Gas Sampling Port:** Yes

**Sterilization:** E.O

**BFE:**99.9999% @ 0.5Microns

**VFE:**99.999% @ 0.027Microns

|         |                             |                             |               |                        |               |
|---------|-----------------------------|-----------------------------|---------------|------------------------|---------------|
| ZF-002P | MOQ:200pcs                  | Single pack in PE-Paper bag | 200pcs/carton | carton size:50x40x32cm | weight:4.1 kg |
| ZF-002P | Single pack in PE-Paper bag | 50pcs/Per Pack              | 250pcs/carton | carton size:50x40x47cm | weight: 5.3kg |