

MINI HMEF FILTER

Mini HME with Bacterial/Viral Filter

● Mini HMEF filter with Port for Adult/Pediatric



RFE:ZF-002

Dead Space:15ml

Size:

H:67mm; w. 44mm

Weight:13.5g

Material:

Housing:

Polypropylene (PP)

Filter: Electrical Cotton

Hydrophobic

Connector:

22M/15F-22F/15M

Tidal volume:100ml-1500ml

Gas Sampling Port: Yes

BFE:99.9999% @ 0.5Microns

Sterilization: E.O

VFE:99.999% @ 0.027Microns

ZF-002P	MOQ:200pcs	Single pack in PE-Paper bag	200pcs/carton	carton size:50x40x32cm	weight:4.1 kg
ZF-002P	Single pack in PE-Paper bag	50pcs/Per Pack	250pcs/carton	carton size:50x40x47cm	weight: 5.3kg