


Adult Anesthesia Circuits (Corrugated)

- Good elasticity, flexibility and air tightness
- Can be sold as kit with anesthesia mask, breathing bag, HMEF, catheter mount, extra limb
- Standard Latex free breathing bags, Latex optional
- ISO standard interface



Anesthesia Circuit

22M/15F



22F
22F

Corrugated


Length	NO.
1.2m	B0001
1.5m	B0002
1.8m	B0003
2.0m	B0823

Corrugated

Length	NO.
1.2m	B0017
1.5m	B0018
1.8m	B0019
2.0m	B0020

Anesthesia circuit, Breathing bag

22M/15F




22F
22F



Anesthesia Circuit, Anesthesia mask

22M/15F



22F
22F



Corrugated


Length	NO.
1.2m	B0041
1.5m	B0042
1.8m	B0043

Corrugated

Length	NO.
1.2m	B0021
1.5m	B0022
1.8m	B0023

Anesthesia circuit, Breathing bag, HMEF

22M/15F




22F
22F



Anesthesia Circuit, Anesthesia mask, Breathing bag

Remark: Breathing bag connector optional 22M/22M

22M/15F



22F
22F



Corrugated


Length	NO.
1.2m	B0080
1.5m	B0081
1.8m	B0082

Corrugated

Length	NO.
1.2m	B0188
1.5m	B0189
1.8m	B0190

Anesthesia Circuit, Breathing bag HMEF, Anesthesia mask

22M/15F



22F
22F




Pediatric Anesthesia Circuits (Corrugated)

- ISO standard interface ID 15mm
- Can be sold as kit with anesthesia mask, breathing bag, HMEF, catheter mount, extra limb
- Standard Latex free breathing bags, Latex optional




Anesthesia Circuit (Corrugated)




Corrugated	
Length	NO.
1.2m	B0224
1.5m	B0225
1.8m	B0226

Corrugated	
Length	NO.
1.2m	B0236
1.5m	B0237
1.8m	B0238




Anesthesia circuit, Breathing bag

Anesthesia Circuit, Anesthesia mask




Corrugated	
Length	NO.
1.2m	B0251
1.5m	B0252
1.8m	B0253

Corrugated	
Length	NO.
1.2m	B0248
1.5m	B0249
1.8m	B0250




Anesthesia circuit, Breathing bag, HMEF

Anesthesia Circuit, Anesthesia mask, Breathing bag



Corrugated	
Length	NO.
1.2m	B0263
1.5m	B0264
1.8m	B0265

Corrugated	
Length	NO.
1.2m	B0269
1.5m	B0270
1.8m	B0271



Anesthesia Circuit, Breathing bag HMEF, Anesthesia mask