

Suction Liner



- ▶ Simplicity
- ▶ Safety
- ▶ Economical



Suction Liner



Product	CAT. NO.	Specification	Size
Suction Liner	6609110	with filter	1000mL
	6609120		1500mL
	6609130		2000mL
	6609210	without filter	1000mL
	6609220		1500mL
	6609230		2000mL

Product features



- 1.This product is used for collecting biohazardous fluids which produced in ICU, operation rooms, gynecology and other clinical practices.
- 2.Double safety design ensures no liquid is suctioned back to tubing and suction unit, effectively prevents cross infection and protects medical staff safe.
- 3.It's made of environment-friendly polyethylene, safe and non-toxic.
- 4.It's easy to use and can be used individually in wall plate or used in tandem in mobile suction hanger (2-4 units).
- 5.Economical, competitive price to help reduce cost.
- 6.Available in some specifications with or without shutoff filter, 1000ml, 1500ml and 2000ml.

