

LUNG EXCERCISER

MODEL : SP-1200
BRAND : RIBBEL,INDIA

Features:

- ✦ Made from transparent material.
- ✦ 3 chamber pf 600,900 and 1200 cc/sec inhaling volume.
- ✦ Three colour coded balls for easier identification of inhaled capacity.
- ✦ Light weights easy to carry.
- ✦ Patient sticker.
- ✦ For single patient use.

Intended Use:

Lung exerciser is used for the measurement of inspiration and expiration capacity of the patient during lung function test and also for the exercise of the lungs / breathing exercise (inspiration)

