

eBfine GDH
Blood Glucose Monitoring System



5 seconds
measuring time



0.5 μ L
Blood volume

180

Results
Memory

Product Specification

Blood sample type	Capillary and Venous Whole Blood
Blood volume	0.5 μ L
Measuring time	5 seconds
Measuring unit	mg/dL & mmol/L (interchangeable)
Acceptable hematocrit Range	20% ~ 60%
Measuring range	20 ~ 600mg/dL (1.1~33.3 mmol/L)
Memory capacity	180 results with time and date
Time display	24H
Operating temperature range	10 ~ 40 °C
Relative humidity operating range	Below 85%
Meter storage condition	0 ~ 50 °C
Meter storage humidity range	Below 95%
Dimension L×W×H(mm)	85×63×15mm
Power supply	CR2032(3V, lithium battery)
Data output	USB PC interface
Weight	≤50g



VISGENEER INC.
www.visgeneer.com

Tel: +886-3-5181918
E-mail : eB@visgeneer.com
Web: www.eB-monitor.com

No. 335, Sec. 6, Zhonghua Rd.,
Hsinchu City 30094, Taiwan



CE 0123

