

eBketone

Blood Ketone Monitoring System



10 seconds
measuring time



0.5µL Capillary
Blood Volume

180

Results
Memory

Product Specification

Blood sample type	Capillary whole blood
Blood volume	0.5µl
Acceptable hematocrit Range	30% ~ 60%
Measuring range	0.0~8.0mmol/L
Measuring unit	mmol/L
Measuring time	10 seconds
Memory capacity	180 results with time and date
Time display	24H
Operating temperature range	10 ~ 40 °C
Relative humidity operating range	Below 85%
Meter storage condition	0 ~ 50 °C
Meter storage humidity range	Below 95%
Dimension L x W x H (mm)	87×63×23mm
Weight	≤80g
Power supply	Two 1.5Volt AAA Batteries
Temperature display	°C



VISGENEER INC.
www.visgeneer.com

Tel: +886-3-5181918
E-mail : eB@visgeneer.com
Web: www.eB-monitor.com

No. 335, Sec. 6, Zhonghua Rd.,
Hsinchu City 30094, Taiwan



CE 0123

